

Product Name :
Karl Fischer Apparatus

Product Code :
NLAB-SCIENTIFICEQ28011



Description :

Karl Fischer Apparatus

Technical Specification :

Microprocessor K.F. moisture titrator with built-in magnetic stirrer, auto computation of moisture in ppm, %age and mg H₂O storage of date, time, batch no. & final results of 20 titrations. K.F. reagent dispensing by motorised process. Final results are flashed on 8 digit display. External printer attachment facility available.

Auto Karl Fischer Titrimeter (Deluxe Model) having separate control and stirring units with array of 10 LED, buzzer auto zero burette, adjustable timer, dual platinum electrode and other accessories.

Auto Karl Fischer Titrimeter with separate control and stirrer unit auto zero burette, dual platinum electrode & other accessories.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003