


Product Name :
Karman Cannula

Product Code :
NLABMEDICAL-M3250001



Description :

Karman Cannula India, Karman Cannula manufacturer India, suppliers, Medical Karman Cannula India, Hospital Karman Cannula India

Technical Specification :

Smooth catheter surface, Rounded Tip.

Smooth edged eye/eyes near the tip for trauma free insertion.

With universal connector for connecting to suction apparatus.

Available in sizes 4mm - 12mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003