

**Product Name :**  
Keithley Nanovolt Meter

**Product Code :**  
NLAB-ENGINEERINGLB30003



**Description :**

Keithley Nanovolt Meter

**Technical Specification :**

Keithley Nanovolt Meter Lab Equipments Manufacturer, Suppliers & Exporters

**Description:**

- Dual-channel Nano voltmeter, Range : 1nV – 100V
- Make low noise measurements at high speeds, typically just 15nV p-p noise at 1s response time, 40nV – 50nV p-p noise at 60ms
- Averages multiple readings for greater noise reduction

- Built-in thermocouple linearization and cold junction compensation

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003