

**Product Name :**  
Kinematic Viscometer Bath

**Product Code :**  
NLAB-BTMENLAB150017



**Description :**

Kinematic Viscometer Bath

**Technical Specification :**

This bath is designed and developed under the strict supervision of our experienced professionals in accordance with the international quality standards. Our offered bath is thoroughly checked on various parameters to ensure zero-defect. This bath is widely used in national laboratories and pharmaceutical industries.

**Features:**

Intelligent signal sampling for accurate data measuring and parameter reading

Accurate viscometer for improving experimental efficiency

High precision digital display

Ideal for controlling the temperature accuracy

Durable glass bath for optimum heat preservation

Desktop based all-in-one design for easy operations and portability

Electric stirrer for ensuring uniform temperature of bath

Easy cleaning and drying of viscometer tubes

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003