

<b>Product Name :</b> Kinetics of Dissolution of Benzoic Acid	<b>Product Code :</b> NLAB-EDUCATIONTVET16033
<b>Description :</b> Kinetics of Dissolution of Benzoic Acid	
<b>Technical Specification :</b> Kinetics of Dissolution of Benzoic Acid -  The setup consists of a Reactor in which solid cylinder of Benzoic Acid is placed and water is stirred by means of a variable speed mixer. Baffles are also provided.  A Peltizer is provided to make the solid cylinder size 1? diameter and 2? height of Benzoic Acid.  The Reactor is made of Stainless Steel.  Water Supply and Drain.  Electricity Supply: 1 Phase, 220 V AC, 0.5 kW.  To study the Kinetics of dissolution of Benzoic Acid.	

Instruments, Laboratory Glassware, and Chemicals required for analysis.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003