

Product Name :
Kitchen Bain Marie

Product Code :
NLAB-VOCATIONALTVET27057

Description :

Kitchen Bain Marie

Technical Specification :

Automation Grade: Semi-Automatic

Material: Stainless Steel

Dimension: 2200X700X850 mm

Voltage: 220 V

Frequency: 50 Hz

Power: 3000 kW

Features:

Attractive design

Easy to use

Accurate size

Long service life

Less upkeep

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003