

Product Name :
Kjeldahl Distillation Instruments

Product Code :
NLAB-EDUSCHOOLEQ31001



Description :

Kjeldahl Distillation Instruments

Technical Specification :

- Kjeldahl Distillation Unit (Without glass parts):
- is useful for distillation purposes in **300 & 500** ml capacity flasks.
- Provided with condenser rack.
- 3-tests mantle type with 3 energy regulators.
- 6-tests mantle type with 6 energy regulators.
- Useful for estimation of Nitrogen contents in substances by digestion in 300 & 500 ml flasks.
- Provided with Lead fume duct and pair of hangers.
- It works on **220V AC 50Hz** single phase.
- Provided with two sets of mantle type heating platforms-one for digestion and other for distillation.

- Complete with condenser rack, lead fume duct and a pair of hangers, on a mild steel sheet tubular stand.
- Both units supplied with individual energy regulators.
- It works on **220V AC 50Hz.**

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003