

**Product Name :**  
Kjeldahl Dry Bath

**Product Code :**  
NLAB-LABORATORYEM27005



**Description :**

Kjeldahl Dry Bath

**Technical Specification :**

Kjeldahl Dry Bath (Tube Heating):- The equipment has a casted top with room for heating 40 digestion tubes of 75 ml each. Heating elements are speared uniformly so that the entire casted top. And thus the tubes attain the same temperature throughout. Temperature is controlled by energy regulators. The gap between the elements and outer M.S. painted wall is filled with high-grade mineral wool to avoid thermal loss. To work on 220 / 230 volts A.C supply.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003