

Product Name :
Klystron Power Supply

Product Code :
NLAB-ELECTRONICSAB120001



Description :

Klystron Power Supply

Technical Specification :

Klystron Power Supply- Engineering Lab Training Systems.

Klystron power supply, is a state-of-the-art solid-state, regulated power supply for operating low power klystrons. It incorporates a number of proprietary features.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003