

Product Name :
Kohlrusch Conductivity Bridge

Product Code :
NLAB-ELECTRICALAB380002



Description :

Kohlrusch Conductivity Bridge

Technical Specification :

A circular slide wire is calibrated from 0.01 to 10 ohms

A rotary dial furnishes six multipliers of 0.1, 1, 10, 100 & 1000 & 10000.

Range of measurement is 0.01 to 100 K Ohm

Terminals for supply and head phone/detectors are provided on the panel.

A detector key is also provided.

The position of slider is adjusted to the position of minimum sound or position of minimum A.C. signals.

Accuracy of coils $\pm 0.1\%$

Accuracy of measurement $\pm 1\%$ at mid scale & $\pm 2\%$ overall.

Optional Accessories:

Kohlrusch Conductivity Bridge

Induction coil

Conductivity Probe With Beaker

Battery Eliminator 2-12V, 2A

Head Phone/ Galvanometer/ Digital Multimeter

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003