

Product Name :
L-Box Apparatus

Product Code :
NLAB -CONCRLAB120018



Description :

L-Box Apparatus

Technical Specification :

The distance between 12 mm diameter bars can be set between 41 ± 1 mm or 59 ± 1 mm. L Shape box is designed for ease of cleaning the vertical and horizontal hoppers.

The I shape box apparatus is supplied complete with:

Filling Hopper

Base

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003