

**Product Name :**  
Lab Automatic Spectrophotometer

**Product Code :**  
NLAB-LABORATORYEM11014



**Description :**

Lab Automatic Spectrophotometer

**Technical Specification :**

The spectrophotometer has become a routine instrument in modern molecular biology laboratories. It is often used for nucleic acid, protein quantification and quantification of bacterial growth concentration.

Automatic detection, accurate measurement

Built in high-speed CPU. 2G memory, fast response, no need to wait

Integration of tablet computer, user-friendly interface, measurement as long as a touch!

All-in-one, working independently, one device meets all requirements

Ultimate industrial design

Precision work, more stable and reliable, simple operation.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments