

Product Name :
Lab Automatic Spectrophotometer

Product Code :
NLAB-LABORATORYEM11014



Description :

Lab Automatic Spectrophotometer

Technical Specification :

The spectrophotometer has become a routine instrument in modern molecular biology laboratories. It is often used for nucleic acid, protein quantification and quantification of bacterial growth concentration.

Automatic detection, accurate measurement

Built in high-speed CPU. 2G memory, fast response, no need to wait

Integration of tablet computer, user-friendly interface, measurement as long as a touch!

All-in-one, working independently, one device meets all requirements

Ultimate industrial design

Precision work, more stable and reliable, simple operation.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments