

**Product Name :**  
Lab Digital Hot Plate

**Product Code :**  
NLAB-LABORATORYEM11028



**Description :**

Lab Digital Hot Plate

**Technical Specification :**

The Digital hot plate is widely used for drying, drying and other temperature tests of samples. It is an indispensable tool for biology, genetics, medicine and health, environmental protection, biochemical laboratories, analysis rooms, education and scientific research.

LED display shows temperature

The highest temperature is up to 550°C

Separate safety circuit, fixed safety temperature is 580°C.

By connecting a temperature sensor (PT 1000) with an accuracy of  $\pm 0.5^{\circ}\text{C}$ , external temperature control can be performed

Glass ceramic plates provide excellent chemical resistance and the most effective heat transfer

Even if the heating plate is turned off, if the temperature of the working plate is higher than 50°C, the "HOT" warning will flash.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003