


**Product Name :**  
Lab Electronic Analytical Balance

**Product Code :**  
NLAB-LABORATORYEM11005



**Description :**

Lab Electronic Analytical Balance

**Technical Specification :**

An analytical balance is the most important instrument for accurate weighing in the experiment, and it can accurately weigh more than 0.0001g (ie 0.1mg).

Single precision mass sensor (external calibration adopts split sensor).

One button automatic internal calibration function (not available for external calibration).

6-level shockproof filter adjustable function.

LCD large interface backlight display.

Weighing stability time is adjustable.

It has 17 unit conversion functions such as gram, carat, pound, ounces, etc.

The standard configuration is RS232 and USB communication interface to connect peripheral devices.

The current mode is automatically locked when the machine is switched on and off.

Built in date and time adjustable function.

Choose the lower hook weighing device.

Built in temperature display function.

Digital multipoint linear calibration.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003