

Product Name :
Lab Electrophoresis Power

Product Code :
NLAB-LABORATORYEM11013



Description :

Lab Electrophoresis Power

Technical Specification :

Electrophoresis is an instrument that realizes electrophoretic analysis. Generally composed of power supply, electrophoresis tank, detection unit, etc. The so-called electrophoresis refers to the movement of charged particles in an electric field. Different substances move at different speeds in the electric field due to their different charges and molecular weights. Based on this, different substances can be analyzed qualitatively or quantitatively, or a certain mixture can be combined. Sub-analysis or individual component extraction preparation.

working state, it can be adjusted in real time

microcomputer intelligent control

liquid crystal display, displaying voltage, current and timing at the same time

switching power supply output

with memory function

with over-voltage, over-current, overload, load changing, no-load and other alarm protection functions

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003