

**Product Name :**  
Lab High-precision Electronic Balance

**Product Code :**  
NLAB-LABORATORYEM11021



**Description :**

Lab High-precision Electronic Balance

**Technical Specification :**

High-precision electronic balance is a weighing instrument used for accurate measurement of mass (weight).

Single precision quality sensor

Automatic intelligent internal calibration function

Automatic dual range, double adjustable function

Six shock filter adjustable function

Anti-interference, anti-static metal shell

LCD large interface VATN display

Weighing stability time adjustable

With grams, carats, pounds, ounces, and other 17 kinds of units conversion function.

Standard configuration RS232 communication interface to connect peripherals.

Switch machine automatically lock the current mode.

Built-in date, time adjustable function.

Built-in temperature display function.

Under the hook weighing device.

Anti-static metal wind shield backplane

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003