

Product Name :
Lab Multichannel Pipette

Product Code :
NLAB-LABORATORYEM11019



Description :

Lab Multichannel Pipette

Technical Specification :

The P series multichannel Micropipette have the light weight, operated with pipetting

Strength and easily and conveniently.

Rotating dispensing head for optimal pipetting convenience.

Due to the optimized standard design for the tip cone which it can be compatible with a wide range hand of tips.

Large central pipetting button and separate ejection function

True one-handed operation for both right-and left-handers

Completely autoclavable at 121°C(20min)

Volume-change protection

4-Position volume display, always clearly visible externally

Shor stroke of only 12.5mm to reduce the risk of RSI(Repetitive Strain Injury)

Corrosion-resistant piston and ejector Color-coded for easy selection of the right tip

Individual shafts with seals with seals can be easily unscrewed for cleaning or replacing-eliminates expense and long outages.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003