

Product Name :
Lab Nanodrop Spectrophotometer

Product Code :
NLAB-LABORATORYEM11024



Description :

Lab Nanodrop Spectrophotometer

Technical Specification :

- I Micro sample detection: the sample on the base only needs 0.5-2 μ L, saving precious samples.
- I Automatic detection: the automatic detection function can be realized by opening and closing the upper arm.
- I Wide detection range: the detection concentration range of the base is 2-27500ng / μ L
- I Easy to operate: built in 7 inch high resolution color capacitive touch screen, friendly operation interface
- I The best choice for nucleic acid detection: including all wavelengths required for nucleic acid concentration and purity detection.
- I Multiple detection modes: it has two detection modes of micro and cuvette, which can easily deal with various types of samples.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003