

**Product Name :**  
Lab Pharmacy Refrigerator

**Product Code :**  
NLAB-LABORATORYEM11006



**Description :**

Lab Pharmacy Refrigerator

**Technical Specification :**

The pharmacy refrigerator is suitable for storage of biological products, vaccines, drugs, reagents, and more. You can apply the lab / pharmacy fridge in pharmaceutical factories, hospital, centers for disease prevention and control, clinics,

100% CFC-free refrigerant design, no chemicals that destroy the ozone layer, which is conducive to environmental protection

Seven temperature probes can ensure high-precision temperature control with almost no fluctuations, thereby improving safety

An indoor lighting system with dual LED lights ensures high visibility inside the cabinet

Provide test port to facilitate users to test the internal temperature of the cabinet

After connecting the U disk, the temperature data can be automatically stored for more than 10 years

Equipped with a USB export interface, which can automatically store data from the previous month to the current month in PDF format

The large capacity of 315L can store vaccines, drugs, reagents and other laboratories/medical materials

conveniently

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003