

Product Name :
Lab Pipet Fillers

Product Code :
NLAB-LABORATORYEM11020



Description :

Lab Pipet Fillers

Technical Specification :

Compatible with most of the plastic and glass pipettes from 0.1 -100mL

Six speeds of selection for aspiration and dispensing different liquids

Large LCD display showing low battery warning and speed settings

Enables single-handed operation with minimum effort

Light and ergonomic design providing easy usability

High capacity Li-ion battery enable long runtime of operation

Powerful pump fills a 25mL pipette in <5 seconds

0.45µm replaceable hydrophobic filter

Rechargeable during use

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003