

**Product Name :**  
Lab Pipet Fillers

**Product Code :**  
NLAB-LABORATORYEM11020



**Description :**

Lab Pipet Fillers

**Technical Specification :**

Compatible with most of the plastic and glass pipettes from 0.1 -100mL

Six speeds of selection for aspiration and dispensing different liquids

Large LCD display showing low battery warning and speed settings

Enables single-handed operation with minimum effort

Light and ergonomic design providing easy usability

High capacity Li-ion battery enable long runtime of operation

Powerful pump fills a 25mL pipette in <5 seconds

0.45µm replaceable hydrophobic filter

Rechargeable during use

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003