

**Product Name :**  
Lab Single Channel Micropipette

**Product Code :**  
NLAB-LABORATORYEM11010



**Description :**

Lab Single Channel Micropipette

**Technical Specification :**

The Single Channel Micropipette is a pipette that can only transfer one liquid sample at a time. We call it a Single Channel Micropipette.

Large central pipetting button and separate ejection function

True one-handed operation for both right-and left-handers

Completely autoclavable at 121°C(20min)

Volume-change protection

4-Position volume display, always clearly visible externally

Shor stroke of only 12.5mm to reduce the risk of RSI(Repetitive Strain Injury)

Corrosion-resistant piston and ejector

Color-coded for easy selection of the right tip

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003