

Product Name :
Lab Single Stage Vacuum Pump

Product Code :
NLAB-LABORATORYEM11031



Description :

Lab Single Stage Vacuum Pump

Technical Specification :

Ultimate deep vacuum: single stage is 10PA, while dual voltage is 0.5 Pa

Compact design with aluminum housing and easy to carry.

Thermal protector in the motor to guarantee the pump to run steadily.

Both 1/4' and 3/8' SAE flare inlet connections allow for flexibility of connections.

Vacuum pump oil included.

Suitable for all the main voltage design in the market.

Individual design for special customers.

Before use, please take off the exhaust cap (the new oil-gas separator has no exhaust cap), and then check the oil level to ensure that the oil level is not lower than the oil level line. If it is lower than the oil level line, oil should be added in time (it is recommended to apply to Tingwei brand vacuum pump oil provided by our company);

Take off the intake cap and connect the extracted container. The pipe should be short, reliable and free of leakage;

Plug in the power plug and turn on the switch to use;

After using, pull off the power plug in time, remove the connecting pipe, and cover the air inlet cap and exhaust cap tightly.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003