

**Product Name :**  
Lab Steam Autoclave

**Product Code :**  
NLAB-LABORATORYEM11002



**Description :**

Lab Steam Autoclave High Quality Laboratory Instrument with Competitive Price. Lab Autoclave Instrument from India China

**Technical Specification :**

The steam autoclave is assembled with a heating system, microcomputer controlled system, over heat and over pressure protection system, which are reliable for the sterilizing effects. convenient of operation and energy conservation. They are ideal equipment for scientific research institutes and hospital laboratories.

All stainless steel structures.

Quick-opening door structure of hand wheel type.

Door safety lock system.

Automatically controlled by computer.

LCD display working status, touch type key.

Automatic protection from overheating and overvoltage.

Steam internal circulation system: no steam emission and the environment to be disinfected should be clean and dry.

Self-expanding seals.

Safety protection against water shortage

Has a perfect drying function.

After disinfection, it will automatically shut down and beep.

With two stainless steel sterilization baskets.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003