

Product Name :
Laboratory Autoclave Semi Autoclave

Product Code :
NLAB-LABORATORYEM47014



Description :

Laboratory Autoclave Semi Autoclave

Technical Specification :

Laboratory Autoclave Semi Autoclave -

Features:-

Fully SS 304 Construction

No Hinge or shaking lids

Reduced loading height

Unique pedal free design

µp based PID controller with dual display for Temp. and Time

Specification:-

Capacity
350 x 500mm
Dimension
SS 304
Construction
Power
Consumption

Operating Pressure	
Operating Temperature	
Settling Time	
Voltage	
Capacity Weight	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003