

Product Name :
Laboratory Bench Oven

Product Code :
NLAB-VOCATIONALTVET18003



Description :

Laboratory Bench Oven

Technical Specification :

Capacity: 220 liters.

Temperature Resolution: $\pm 0.10\text{C}$.

Rated power heater elements: 2.4 KW

Temp. range: +50-300°C.

Time Locating Range: 1-9999Min.

Temp. Fluctuation: $\pm 0.75\ 0\text{C}$.

Shelves: 2, Shelf positions; 3

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com