

**Product Name :**  
Naugra Lab Laboratory Mill

**Product Code :**  
NLAB-LABORATORYEM65024



**Description :**

Naugra Lab Laboratory Mill

**Technical Specification :**

Laboratory Mill (Willey type) :-

For preparation of samples for analysis with minimum loss of moisture of volatile constituents.

Weight distribution and balancing of rotating parts has been adjusted to provide greater stability and quieter performance with reduced vibration.

Suitability for wider variety of commercial materials. fitted with 1/4 HP motor, supplied with nickel plated hopper of 50 mm diameter with cover, three delivery sieves

Having tops of 20, 40 and 80 mesh respectively, with plastic receiver, brush, L-N key, cord and plug.

Suitable to operate on 220V, 1 Phase, 50 Cycles AC supply.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003