

Product Name :
Laboratory Pan Mixer

Product Code :
NLAB -CONCRLAB120012



Description :

Laboratory Pan Mixer

Technical Specification :

The shutter of the mixer is designed to open 120 degrees for easy access to the drum and keeps the beater in a certain height during mixing. There is a small observation window on the shutter which enables the user to monitor the process.

Laboratory pan mixer

The mixing drum can be tilted to empty the mixture without hassle and it can be removed for easy cleaning on completion of the mixing operation. It is equipped with rubber wheels which provides high mobility.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003