

**Product Name :**  
Laboratory Pasteuriser

**Product Code :**  
NLAB-EDUCATIONTVET16006

**Description :**

Laboratory Pasteuriser

**Technical Specification :**

Laboratory Pasteuriser -

The unit may also be used for project work, in particular, to measure the effects on product quality of different process conditions.

The Laboratory Pasteuriser is a benchtop unit designed to provide practical training for the operation of a large modern pasteurizing plant.

**Features:-**

Miniature plate heat exchanger

Control console incorporating comprehensive instrumentation.

Variable product temperature and holding time

Choice of holding tubes

Fully self-contained design with integral hot water system

**Applications:-**

Arresting of yeast activity

Quality control

Destruction of spoilage organisms

Arresting of enzyme activity

HTST pasteurization of milk Destruction of spoilage organisms

**Benefits:-**

Tests can be carried out quickly and easily

Benchtop operation

Simulates the conditions used in a production plant

Small quantities of liquid product can be processed

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003