

<b>Product Name :</b> Laboratory Pasteuriser	<b>Product Code :</b> NLAB-EDUCATIONTVET16006
<b>Description :</b> Laboratory Pasteuriser	
<b>Technical Specification :</b> Laboratory Pasteuriser -  The unit may also be used for project work, in particular, to measure the effects on product quality of different process conditions.  The Laboratory Pasteuriser is a benchtop unit designed to provide practical training for the operation of a large modern pasteurizing plant.	
<b>Features:-</b>  Miniature plate heat exchanger	

Control console incorporating comprehensive instrumentation.

Variable product temperature and holding time

Choice of holding tubes

Fully self-contained design with integral hot water system

**Applications:-**

Arresting of yeast activity

Quality control

Destruction of spoilage organisms

Arresting of enzyme activity

HTST pasteurization of milk Destruction of spoilage organisms

**Benefits:-**

Tests can be carried out quickly and easily

Benchtop operation

Simulates the conditions used in a production plant

Small quantities of liquid product can be processed

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003