

Product Name :
Laboratory Plasticwares

Product Code :
NLAB-EDUSCHOOLEQ168001

Description :

Laboratory Plasticwares Educational

Technical Specification :

These unbreakable Scintillation vials, moulded in Polyethylene, are fitted with leak proof screw caps.

These Scintillation vials can be used with standard counting equipment for liquid scintillation system and are excellent for low activity count.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003