

Product Name : Laboratory Polarising Microscope	Product Code : NLAB-EDUSCHOOLEQ1730132
Description :	
Laboratory Polarising Microscope	
Technical Specification :	
<p>The microscope stand has a quadruple ball bearing nose-piece, co-axial coarse and fine motions moving on ball bearing guide. Sub-stage polariser on strain-free Abbe condenser Numerical Aperture of 1.25.</p> <p>Polarising Microscopes are standard polarising microscopes equipped with strain-free plan achromatic objectives for transmitted light observation. Built-in halogen illumination system and high quality coated optics make the equipment versatile and high quality instrument.</p> <p>Polarising intermediate has attachment magnification of 1.2X and has built-in Bertrand lens unit on a rotating turret for colonoscopy, a swing in / swing out analyser and DIN standard compensatory slot accommodating 1st order red sensitive tint plate and $\frac{1}{4}$ wave length radiation plate.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003