

Product Name :
Laboratory Polarising Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730132

Description :

Laboratory Polarising Microscope

Technical Specification :

The microscope stand has a quadruple ball bearing nose-piece, co-axial coarse and fine motions moving on ball bearing guide. Sub-stage polariser on strain-free Abbe condenser Numerical Aperture of 1.25.

Polarising Microscopes are standard polarising microscopes equipped with strain-free plan achromatic objectives for transmitted light observation. Built-in halogen illumination system and high quality coated optics make the equipment versatile and high quality instrument.

Polarising intermediate has attachment magnification of 1.2X and has built-in Bertrand lens unit on a rotating turret for colonoscopy, a swing in / swing out analyser and DIN standard compensatory slot accommodating 1st order red sensitive tint plate and $\frac{1}{4}$ wave length radiation plate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003