

**Product Name :**  
Laboratory Stirrer

**Product Code :**  
NLAB-EDUSCHOOLEQ45005



**Description :**

Laboratory Stirrer

**Technical Specification :**

Heavy & stable cast iron base duly painted with mild steel chrome plated upright rod.

AC/DC motor with stainless steel stirring shaft and blades, variable speed control unit provided.

1/20 HP motor provided for smooth long-time stirring.

Heating of stirrer and stirring space adjustable.

To work on 220 V AC 50 HZ single phase.

Automatic 0-60 minutes timer.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

