

Product Name :
Laminar Air Flow Cabinet Stainless Steel

Product Code :
NLABMEDICAL-M3160007

Description :

Laminar Air Flow Cabinet Stainless Steel India, Laminar Air Flow Cabinet Stainless Steel manufacturer India, suppliers, Medical Laminar Air Flow Cabinet Stainless Steel India, Hospital Laminar Air Flow Cabinet Stainless Steel India

Technical Specification :

Laminar Air Flow Cabinet (Stainless Steel)

Table Top model.

Satin Finish S.S.Table Top(in lieu of formica Top).

Transparent front door.

Static Pressure. Inclined Manometer.

Built in U.V.Germicidal light.

Cock for gas, air or vacuum line.

Cold rolled MS Sheet body duly powder coated.

Stainless steel body duly polished.

Particle dust count validation at site after installation.

HEPA Filter.

Pre Filter.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments