

Product Name :
Laminar Air Flow Horizontal With Stand

Product Code :
NLAB-LABORATORYEM26002



Description :

Laminar Air Flow Horizontal With Stand

Technical Specification :

Standard Size (Working Area): 2' x 2', 3' x 2', 4' x 2', 6' x 2', 8' x 2'

Air Velocity: 90 ± 20 cfm (Average)

Filtration Elements: Main Filters : HEPA Filter Prefilter: Washable non-woven polyester fibers.

Filter Efficiency: Main Filters : 99.97 at 0.3 micron Prefilter: Arrestance 90%, at 10 microns

Work: Stainless Steel

Noise Level: < 65 dBA

Light Intensity: 800-1000 Lux

Power Supply: 220V Single Phase, 50 Hz.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003