

**Product Name :**  
Large Tool Makers Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730127

**Description :**

Large Tool Makers Microscope

**Technical Specification :**

Objective : 2x achromatic.

Illumination : Sub-stage lamp emits transmitted light from a bottom source equipped with collimated green filter halogen light. Oblique illumination through LED ring.

Stage : 300mm X 200mm thread less stage having built-in liner scales, least count 0.001mm, to provide accurate and smooth X-Y table travel 100x 50mm. Glass plate size 200 X 100mm.

Eyepiece : WF 15x with cross reticle

Base : Large and heavy base provides extra overall rigidity to the instrument.

Magnification: 30x standard.

Excellent for Precise Measurements of Angles, Lengths, Diameters, and point to point distances on the specimen.

Observation Tube : Monocular, inclined at 30° coupled with graduated eyepiece protractor.

Compact design, especially engineered for all who require precision and accuracy for measuring a variety of metal, plastic, electronic parts and automotive components.

Eyepiece Protector : Graduated in 360° with adjustable vernier reading to 6 Min. is coupled with monocular ocular tube for smooth angle measurement.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments