

Product Name :
Lateral Pressure Assembly

Product Code :
NLAB-TESTINGLABM109057



Description :

Lateral Pressure Assembly

Technical Specification :

Lateral Pressure Assembly 10kg/cm²

It is designed for application of cell pressure upto 10kg/cm on a triaxial specimen in the triaxial cell. The unit is mostly used in routine tests which last for about 10 to 15 minutes wherein the variations of cell pressure can be tolerated.

Specification:

The apparatus consists of a pressure chamber which has a flanged top 2 cap fitted with a 10 kg/cm² . Pressure gauge and a valve for pumping in air. Water inlet and drain cock are fitted to the chamber. 2 The foot pump supplied easily develops a pressure of 10 kg/cm² Complete with connecting pressure hose.

Lateral Pressure Assembly 20kg/cm²

Same as above but suitable to develop pressure upto 20kg/cm fitted with suitable pressure gauge.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003