

**Product Name :**  
Lauderometer Testing Machines

**Product Code :**  
NLAB-TESTINGLABM75001



**Description :**

Lauderometer Testing Machines

**Technical Specification :**

The Lauder-Ometer is a laboratory instrument used for conducting accelerated laundering or dyeing tests. It is applicable for performing washing, dry cleaning, dyeing, colorfastness to fulling, detergency, and other laboratory tests. It is equipped with a four-sided, stainless steel rotor, which holds five metal containers on each side. The containers in the rotor are simultaneously exposed to controlled conditions of test solution, temperature, and mechanical agitation. Each container is sealed and allows a reliable experiment with no interaction from other test specimens

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003