

Product Name : Launderometer	Product Code : NLAB-TESTINGLABM59003
--	--

--

Description :

Launderometer

Technical Specification :

Specifications

Rotation speed 40 +/- 2rpm

AATCC canister 1200ml, 6 pcs

ISO canister 550ml, 12 pcs

Max 98 oC, adjustable

Running time adjustable

Steel Balls 200 pcs

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003