

**Product Name :**  
LCD Screen Digital Handheld Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730105

**Description :**

LCD Screen Digital Handheld Microscope

**Technical Specification :**

Focus Range: Manual (10mm – 50mm).

Battery: Rechargeable lithium battery.

Frame Rate: Max 30f/s under 600 Lus Brightness.

Magnification Ratio: 10x to 300x, up to 1200x by digital.

Bundle Software: Portable Capture Pro with measurement.

Still/Video Capture Resolution: 12M, 9M, 5M, 3M, 1.3M, VGA.

Measurement: By software (PortableCapture with OSD languages).

Card Slot: Supports MicroSD Card slot up to 32GB SDHC.

Light Source: 8 LED.

TV Output: Available.

Image Sensor: 5 mega pixels.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003