


Product Name :
Le Chatelier Kettle

Product Code :
NLAB-EDUCATIONTVET15046



Description :

Le Chatelier Kettle

Technical Specification :

Le Chatelier bath for cement testing

Inner chamber in stainless steel

For 12 Le Chatelier molds (not included)

Bath boiling point reached in 30 minutes

Voltage: 220 V / 50 Hz

Power: 1800 W

Dimensions: 405x265x205 mm

Weight: 7 kg

Le Chatelier needle mold according to EN196, numbered

Glass plate to cover the mold (1 pair)

100g mass

Apparatus for checking the elasticity of the mold

Stitching rod diam. 17 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003