

Product Name :
Leak Proofness Tester Machines

Product Code :
NLAB-TESTINGLABM62001



Description :

Leak Proofness Tester Machines

Technical Specification :

Break Proofness Tester - The assemblage of waterproof defensive footwear or gumboots might be secured against entrance of water experienced when the wearer strolls through water. The capacity of waterproof defensive footwear to prevent water from entering them is dictated by shutting their best opening with a watertight metallic or elastic end plug and drenching it up to its best edge submerged. Compacted air at a predefined weight is presented in the footwear and any spillage of air from it into water is inspected outwardly for a steady stream of air bubbles. The footwear under test is kept vertical inside the water shower utilized for its inundation by setting various huge measurement steel balls in it.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003