

Product Name :
LED Characteristics Trainer Digital

Product Code :
NLAB-ELECTRICALAB210051



Description :

LED Characteristics Trainer Digital

Technical Specification :

Digital Meters:

Ammeter 20mA DC.

Voltmeter 20V DC.

Power supplies:

DC Supply IC Regulated 0-5V DC, 150mA.

Components are mounted on the panels are:

LEDS:

RED

YELLOW

BLUE

GREEN

Voltage Control through Potentiometer.

Series Resistor 470E

Salient Features:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003