

Product Name :
LED Cordless Beginner Monocular Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730119

Description :

LED Cordless Beginner Monocular Microscope

Technical Specification :

This beginner microscope offers great value and performance. This microscope can be fully charged in 8 hours, giving it up to 50 hours of continuous use. Excellent for young children who are studying biology and those just starting out.

Monocular head, 45° inclined and 360° rotatable.

Locked-on 10X wide-field eyepiece with pointer.

Coarse focus control.

Rechargeable LED cordless illumination with on/off switch.

4X, 10X, and 40XR, color-coded achromatic objectives are parfocaled and parcentered.

It features an all metal frame, making it rugged and durable, objectives that are achromatic, color-coded and parcentered for a sharp, clear image, and LED illumination that provides cool-to-touch, crisp white lighting.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003