

**Product Name :**  
LED Kohler Illumination Monocular Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730113

**Description :**

LED Kohler Illumination Monocular Microscope

**Technical Specification :**

Eyepiece with a calibrated pointer ensures that definite points on the slide can be easily designated, identified and measured.

Double layer, built-in mechanical stage is low-profile and ball-bearing, with a mounted specimen holder, ensuring easy mobility and precise control of specimen under observation.

Corded LED illumination with light intensity control and adjustable condenser diaphragm for Kohler illumination.

40X to 1000X magnification.

Monocular head is 45° inclined and rotates 360°.

Rack and pinion adjustable 1.25 N.A. Abbe Condenser, and turn-wheel opening iris diaphragm with plastic shield protector, for optimal resolution and image contrasting.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003