

Product Name : LED Kohler Illumination Monocular Microscope	Product Code : NLAB-EDUSCHOOLEQ1730113
Description :	
LED Kohler Illumination Monocular Microscope	
Technical Specification :	
<p>Eyepiece with a calibrated pointer ensures that definite points on the slide can be easily designated, identified and measured.</p> <p>Double layer, built-in mechanical stage is low-profile and ball-bearing, with a mounted specimen holder, ensuring easy mobility and precise control of specimen under observation.</p> <p>Corded LED illumination with light intensity control and adjustable condenser diaphragm for Kohler illumination.</p> <p>40X to 1000X magnification.</p> <p>Monocular head is 45° inclined and rotates 360°.</p> <p>Rack and pinion adjustable 1.25 N.A. Abbe Condenser, and turn-wheel opening iris diaphragm with plastic shield protector, for optimal resolution and image contrasting.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003