

Product Name : LED Phototherapy Stand	Product Code : NLABMEDICAL-M3440002
Description :	
LED Phototherapy Stand	
Technical Specification :	
LED's emitting blue light in 450-470 nm.	
Spectrum matching peak absorption wavelength (458 nm.) at which bilirubin is broken down	
Do not emit light in infrared radiation range reduction the risk of fluid loss	
Long life > 20,000 hrs	
Negligible Maintenance cost	
Intensity: > 25 μ W / cm ² / nm	
Heavy duty mobile stand with hight adjustment facility & easy swivelling box when stand is to be used in conjunction with radiant warmer or infant trolley	
No need to Cover eyes of the baby	
Mounted on casters for easy mobility	
M. S. Fabrication with powder coated paint	
Time Totalizer (Optional)	
Operation Voltage: 100-250V AC, 50Hz	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

