

**Product Name :**  
LED Phototherapy Stand

**Product Code :**  
NLABMEDICAL-M3440002

**Description :**

LED Phototherapy Stand

**Technical Specification :**

LED's emitting blue light in 450-470 nm.

Spectrum matching peak absorption wavelength (458 nm.) at which bilirubin is broken down

Do not emit light in infrared radiation range reduction the risk of fluid loss

Long life > 20,000 hrs

Negligible Maintenance cost

Intensity: > 25  $\mu$ W / cm<sup>2</sup> / nm

Heavy duty mobile stand with height adjustment facility & easy swivelling box when stand is to be used in conjunction with radiant warmer or infant trolley

No need to Cover eyes of the baby

Mounted on casters for easy mobility

M. S. Fabrication with powder coated paint

Time Totalizer (Optional)

Operation Voltage: 100-250V AC, 50Hz

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

